

Scan, Verify, Protect - Safe4Use Simplifies Surgical Safety.

Introducing Safe4Use - The Smarter Way to Safeguard Surgical Care

Safe4Use, the latest innovation from Medical Intelligence Systems, is transforming how healthcare professionals verify the safety of medical devices before use in surgical procedures.

With a simple scan of a QR code using any smart device, clinicians can instantly confirm whether a tray, instrument, or device is **within its validated expiry window** and **safe for patient use**. No more guesswork. No more missed checks. And crucially – no more preventable clinical episodes caused by expired equipment.

Key Benefits:

- **Instant Verification:** Scan and receive real-time confirmation of device status.
- **Patient-First Safety:** Ensures every item used in theatre meets sterility and expiry standards.
- **Streamlined Workflow:** Reduces manual checks and reprocessing delays.
- **Audit-Ready:** Creates a digital trail for compliance and traceability.

How It Works:

1. **Scan the QR code** on the medical device or tray set.
2. **Safe4Use instantly checks** against expiry data and usage logs.
3. **Clear status alert:** Green for safe, red for expired – no ambiguity.
4. **Proceed with confidence**, knowing every item meets clinical safety standards.

Whether you're in a busy NHS theatre or a private surgical suite, **Safe4Use** empowers teams to uphold the highest standards of care – with speed, clarity, and compassion.



Join the  Revolution!